

Philips Psoriasis Treatment BlueControl



PSK0202

Effective, wearable relief from plaque psoriasis

# With UV-free blue LED-light

With Philips BlueControl, you decide when and where to treat your plaques. This gentle, UV-free therapy gives you the clinically proven benefits of blue LED light without interrupting your normal daily routine. Try it for yourself at home.\*

# First blue light therapy to treat plaque psoriasis

• World's first blue light therapy to treat plaque psoriasis

# Gentle UV-free therapy without chemicals

- Gentle, UV-free light therapy
- No chemical substances are used

# Convenient, wearable medical device

- · Convenient, wearable relief
- Alleviates plaque psoriais while you get on with your life

## Significant improvement of plaque symptoms

- Blue light induces natural processes
- Clinically proven therapy



# Highlights

# Innovative blue LED light



Philips BlueControl is the world's first wearable phototherapy for mild to moderate plaque psoriasis. Thanks to the innovative design of the small battery-driven device, you can benefit from clinically proven blue LED light while following your normal daily routines.

# LEDs without UV-radiation



Blue LED light belongs to the visible light range and is therefore free of the risks of ultraviolet (UV) radiation. The 40 integrated blue LEDs in Philips BlueControl are long-lasting and energy-efficient sources of light. The LED settings have been optimized for the treatment of plaque psoriasis through in-depth research carried out across several years.

## Therapy without chemicals



Studies confirm that blue LED light, as used in Philips BlueControl, is not toxic for the skin.\* The blue light treatment induces natural processes in your skin helping to alleviate plaque symptoms. The therapy is gentle as it does not make use of any chemical substances.

# Choose when and where to treat



You decide where and when to treat your plaques with the cordless device. An accompanying textile fixation strap attaches Philips BlueControl securely to your limbs, and is adjustable for a comfortable fit during treatment and light activities.

#### Easy to conduct treatment



Experience for yourself how easy it is to get great results without disrupting your daily life. Just attach the device to your arm or leg, press the button ... and carry on with your daily activities. Philips BlueControl shuts off automatically after 30 minutes' treatment per plaque – no timer required!

#### Helps skin to renew properly



Blue light improves plaque symptoms by inducing natural processes in the skin. It helps skin to renew itself properly by slowing the accelerated production of skin cells associated with plaque psoriasis.\* With its antiinflammatory properties, blue light controls unregulated inflammation and can alleviate the symptoms of affected skin.\*

#### **Proven improvement of plaques**

Original Pape					Original Pape	iginal Paper	
Dermotology Commission	111		ORGAN A	Demolology	in an and	All North States 1, 21 North States and R. 21 North States and R. 21 North States and R. 21	
Prospective Randomized Study on the Efficacy of Blue Light in the Treatment Psoriasis Vulgaris		Blue-Light Irradiation Regulates Proliferation and Differentiation in Human Skin Cells		Prospective Randomized Long-Term Study on the Efficacy and Safety of UV-Free Blue Light for Treating Mild Psociasis Vulgaris			
Antonia Warnatahi Spaarene Hoff Lench Hars J. Kierk Revena opolisiber Inganisas at Januaring and Angeleg. Smanly Hopks, Millinghan (Sanat), Jackes, Jan		Leng Latercard', Matthian Boorf' and Fiscaria Scale Scale Scale of the Scale Scal		Souphanie Malffer, Jong Ladonann <sup>6</sup> . Marthiae Bonh <sup>6</sup> . Hars A. Mark <sup>6</sup> . Vannes visie Andere Vannes visien Andere and Antoniae, source, maakki Mith Autor source, aar <sup>100</sup> Aprilant, Invester Mahange, Kana Januer			
y Verb La (g): - Rolasis valgets: - Propicities and initial into high Tarlasis valgets: - Propicities and initial into high Tarlasis	Reference Periodical price (PT) & a order characterised by on the	altibutels to differentiation reductors is done by an demonstrate the literalized instalation up to elliver phe- proteins, which is known to another dimensionless or mechanism for which these light may be alticulate in a proteinsterior data to free reductors of differentiation. No strategies and the literal set of differentiation. No strategies and the literal set of the analysis of the proteinsterior data for the reductors of differentiation. No strategies and the set of the analysis of additionation of an analysis.	nonces of differentiation conducts. Spectrosoft a conductably generates minit could infill how or date code. Set data provide estimates for a of meeting, happenedificative data conditions for a meeting.	Section 1		Internet	
	ting doub 5.121 to scheriose light doubed door dependent tion in off member top cannot	<text><text><text><text><text></text></text></text></text></text>	An extra contract on the second secon	Chicken Charles Analysis of the Chicken Charles and the Chicken Chic	In Operandary and the second s	Answering states of the latest and statest through the statest of the statest	
COLUMN TO AND ADDRESS OF	Sama and San Area Sama and Sama and Sama and Sama and Sama and Sama and Sama Sama and Sama and Sama and Sama Sama and Sama and Sama and Sama Sama and Sama and Sama and Sama and Sama and Sama Sama and Sama and Sama Sama and Sama and Sama Sama and Sama and Sama Sama and Sama and Sama Sama and Sama and Sama Sama and Sama and S	discusses of phoneing but as senses of an out- trin should be an out-time bits sensed if any list is call be basis to basispine binantig	able scribply pairs to beings who as how allowed to be the scribbly as off or behave 10 of bit as he is due to as off or integers 10 of bit as he is due to as off or	1. Starting of		United Rept. 410 (1999) Standards (1999) State Trianel Lange and State Rept. 4	

Clinical studies have confirmed the positive effects of blue LED light on psoriasis plaque symptoms such as redness, scaling and thickness.\* Over 12 weeks of treatment plaque symptoms improved significantly for patients in the clinical study. In some patients, a complete clearance of the plaque was observed. During the treatment-free follow-up phase over 1 month, plaque symptoms improved further to a small extent, rather than worsening. 83% of patients rated the usability and comfort of the device as ideal or excellent.\*

# **Specifications**

## **Accessories**

- USB cable: Recharging
- · Charging adapter: Recharging
- Fixation strap: Secure fixation

# **Technical specifications**

- Frequency: 50 60 Hz
- Type of lamps: LED
- Voltage: 100 240 V
- Power consumption: 15 W

# Safety and Regulations

- CE 0344 Certified
- UV-free: No UV or near UV radiation
- Optical output classification: Risk Group 1 (EN 60601-2-57)
- Power density: 200 mW/cm<sup>2</sup> (peak), 40 mW/cm<sup>2</sup> (average)

- Max. daily dose: 90 J/cm<sup>2</sup> per plaque
- Medical device class lia
- MDD 93/42/EEC
- Logistic data
- Country of Origin: The Netherlands

#### **Clinical trials**

• Efficacy shown in trials: Clinically proven

#### **Dimensions and weight**

- Product dimension: W/H/L: 77 x 23 x 125 mm
- Product weight: 165 grams

#### **Portable**

- Rechargeable battery
- Small, compact design • Power adapters (EU and UK)

DHIIDQ

Issue date 2015-07-30

© 2015 Koninklijke Philips N.V. All Rights reserved.

Version: 2.0.1

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

12 NC: 8670 001 32907

www.philips.com

\* Philips BlueControl is a Medical device class IIa (CE0344) for the \* Philips BlueControl is a Medical device class lia (CE0344) for the treatment of mild to moderate psoriasis vulgaris. Before altering or stopping any prescribed medication or using the device, patients should consult with their physician and carefully read the instructions for use. March 2015.
\* Blue light as used in Philips BlueControl is not toxic for the skin -Awakowicz P et al. Contributions to Plasma Physics. 2009; 49 (9):

- 641 647.
- \* Irradiation with blue light reduces the production of skin cells -Irradiation with blue light reduces the production or skin cens -Liebmann J, Born M, Kolb-Bachofen MV. Journal of Investigative Dermatology. 2010; 130: 259 - 269.
  Blue light induces anti-inflammatory processes - Fischer M et al. Experimental Dermatology. 2013; 22: 554 - 563.
  The positive effects of blue LED light were confirmed in clinical trials - Weinstabl A et al. Dermatology. 2011; 223 (3): 251 - 9.

- \* Over 12 weeks of treatment the plaque symptoms reduced for 84% of the patients. For some patients the blue light even cleared their treated plaque completely Pfaff S et al. Dermatology 2015; 231: 24 - 34